







## Sensor & Cable Compatibility

Technology ▶	Masimo SET® or Masimo rainbow SET™		
	Cable ▶	M-LNC/rainbow® RC	RD SET™/RD rainbow SET™
Connector ▶	DB9	M15	RD15
			
Sensor			
LNCS® Disposables	Y	N	N
LNCS Reusables	Y	N	N
M-LNCS™ Disposables	N	Y	N
M-LNCS Reusables	N	Y	N
RD SET Disposables	N	N	Y
RD SET Reusables	N	N	Y

# Minimum Software Compatibility Requirements\*

Parameters ▶	SET® Parameters					
	Cable ▶		Connector ▶		Sensor	
	LNC	Red LNC	M-LNC	rainbow® RC	RD SET	RD SET MD20
	DB9		M15		RD15	
						
LNCS Disposables	No Min	7.0				
LNCS Reusables	No Min	7.0				
M-LNCS Disposables			No Min	7.1		
M-LNCS Reusables			No Min	7.1		
RD SET Disposables					No Min	7.1
RD SET Reusables					No Min	7.1
Trauma and Newborn	4.1	7.0	4.1	7.1	4.1	7.1

Compatible
  Not compatible

\* Masimo technology board software, sensor, and sensor cable combinations have been validated and tested for compliance per ISO 80601-2-61.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

**Masimo Americas**  
 tel 1-877-452-7455  
 info-america@masimo.com

**Masimo International**  
 tel +41-32-720-1111  
 info-international@masimo.com

