

rainbow® Sensor & Cable Compatibility

Masimo Technology ▶		rainbow SET™		SET®		
Patient Cable ▶		RD rainbow SET™	rainbow® RC	RD SET™	M-LNC	Red LNC/LNC
rainbow® Sensor ▼	Connector ▶	RD15	M15	RD15	M15	DB9
RD rainbow SET-2 Disposables		Y	N	Y	N	N
RD rainbow Lite SET™-1Disposables		Y	N	Y	N	N
rainbow® Disposables R		N	Y	N	Y	N
rainbow® Disposables R1		N	Y	N	Y	N
rainbow ReSposable® R2		N	Y	N	Y	N
rainbow® SpHb® Spot-check Reusables		N	Y	N	N	N
rainbow® SpCO® Reusables		N	Y	N	Y	N

Minimum Software Compatibility Requirement[†]

Patient Cable ▶		RD rainbow SET™					RD SET™		
rainbow® Sensor ▼	Parameter ▶	SpHb®	SpMet®	SpCO®	ORi™	PVi®	SET®	PVi	SET®
RD rainbow SET-2 Disposables		7.4	7.1		7.C	4.8	No Min	4.8	No Min
RD rainbow Lite SET™-1 Disposables					7.C	4.8	No Min	4.8	No Min

Patient Cable ▶		rainbow® RC					M-LNC		
rainbow® Sensor ▼	Parameter ▶	SpHb	SpMet	SpCO	ORi	PVi	SET®	PVi	SET®
rainbow® Disposables R			7.1	7.1		4.8	No Min	4.8	No Min
rainbow® Disposables R1		7.4	7.1		7.C	4.8	No Min	4.8	No Min
rainbow ReSpensible® R2		7.4	7.1		7.C	4.8	No Min	4.8	No Min
rainbow® SpHb® Spot-check Reusables*		7.4	7.4	7.4**		7.4	7.4		
rainbow® SpCO® Reusables			7.1	7.1		4.8	No Min	4.8	No Min

[†] Masimo technology board software, sensor, and sensor cable combinations have been validated and tested for compliance per ISO 80601-2-61.

* Device must have SpHb parameter enabled.

** SpHb and SpCO measurements can be obtained simultaneously using the rainbow® Super DCI-mini sensor.

ORi, rainbow® Super DCI-mini sensor, and RD rainbow Lite SET-1 sensors have obtained CE Marking. Not available in the U.S. or Canada. RD rainbow SET Sensors are not licensed for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo Americas
tel 1-877-452-7455
info-america@masimo.com

Masimo International
tel +41-32-720-1111
info-international@masimo.com

