




# NomoLine® Device Compatibility Chart

The Masimo NomoLine family of sampling lines is compatible with many existing capnography and gas monitoring devices.

Masimo Sampling Line	Compatible Devices with Examples
 <p>NomoLine</p>	<p><b>Masimo Devices such as:</b></p> <ul style="list-style-type: none"> <li>&gt; Root® with ISA™</li> <li>&gt; Rad-97™ with NomoLine Capnography</li> </ul>
 <p>NomoLine-O™</p>	<p><b>Medtronic® (Oridion) Microstream Devices</b></p> <ul style="list-style-type: none"> <li>&gt; Capnostream™ 35</li> <li>&gt; Microstream™ EtCO2 Modules</li> </ul>
 <p>NomoLine with Luer</p>	<p><b>Devices with Water Traps</b></p> <ul style="list-style-type: none"> <li>&gt; GE® with Datex Ohmeda</li> <li>&gt; Dräger®</li> </ul>

The Medtronic®, GE®, and Dräger® trademarks are property of their respective owners and are used for identification purposes only. Use of these names, trademarks, or brands does not imply any affiliation with or endorsement.

For professional use. See instructions for use for full prescribing information including indications, contraindications, warnings and precautions. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

**Masimo Americas**  
tel 1-877-452-7455  
info-america@masimo.com

**Masimo International**  
tel +41-32-720-1111  
info-international@masimo.com

© 2019 Masimo. All Rights Reserved.

